

CFAP69 Polyclonal Antibody

Catalog No: #30352



Package Size: #30352-1 50ul #30352-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CFAP69 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant protein of human CFAP69.
Accession No.	Swiss-Prot#:A5D8W1NCBI Gene ID:79846
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Background

Published Papers

el at., Peptides from the croceine croaker (*Larimichthys crocea*) swim bladder attenuate busulfan-induced oligoasthenospermia in mice.

In *Pharm Biol* on 2022 Dec by Xi Zha, Mengmeng Sang, et al..PMID:35148224, , (2022)

[PMID:35148224](https://pubmed.ncbi.nlm.nih.gov/35148224/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.