

SAMHD1 (Phospho-Thr592) Antibody

Catalog No: #12805



Package Size: #12805-1 50ul #12805-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | SAMHD1 (Phospho-Thr592) Antibody |
| Purification | Affinity purification using immunogen. |
| Applications | WB |
| Species Reactivity | Hu |
| Specificity | SAMHD1 (Phospho-Thr592) Antibody detects endogenous levels of SAMHD1 only when phosphorylated at Thr592 |
| Immunogen Type | Peptide |
| Immunogen Description | Synthesized peptide derived from human SAMHD1 around the phosphorylation site of Thr592. |
| Target Name | SAMHD1 (Phospho-Thr592) |
| Modification | Phospho |
| Other Names | SAMHD1, AGS5, CHBL2, DCIP, HDDC1, Monocyte protein 5, MOP-5, MOP5, SBBI88, Mg11, SAM domain and HD domain 1 |
| Accession No. | Swiss-Prot#: Q9Y3Z3NCBI Gene ID: 25939 |
| Target Species | human |
| Calculated MW | 72KD |
| Concentration | 1.0mg/mL |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at -20°C |

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.