

β -Catenin (Phospho-Ser33/37/Thr41) Antibody

Catalog No: #12806



Package Size: #12806-1 50ul #12806-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	β -Catenin (Phospho-Ser33/37/Thr41) Antibody
Purification	Affinity purification using immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	β -Catenin (Phospho-Ser33/37/Thr41) Antibody detects endogenous levels of β -Catenin only when phosphorylated at Ser33/37/Thr41
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human β -Catenin around the phosphorylation site of Ser33/37/Thr41.
Target Name	β -Catenin (Phospho-Ser33/37/Thr41)
Modification	Phospho
Other Names	CTNNB1, Beta-catenin, Catenin beta-1, CTNNB, Beta Catenin
Accession No.	Swiss-Prot#: P35222NCBI Gene ID: 1499
Target Species	human
Calculated MW	92KD
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.