

# JNK1/2 Polyclonal Antibody

Catalog No: #30394



Package Size: #30394-1 50ul #30394-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	JNK1/2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant protein of human JNK1/2.
Accession No.	Swiss-Prot#:P45983/P45984NCBI Gene ID:5599/5601
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.