

VIPAS39 Polyclonal Antibody

Catalog No: #30460

Package Size: #30460-1 50ul #30460-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	VIPAS39 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human VIPAS39 (NP_071350.2).
Other Names	VIPAS39; C14orf133; SPE-39; SPE39; VIPAR; VPS16B; hSPE-39; spermatogenesis-defective protein 39 homolog
Accession No.	Swiss-Prot#:Q9H9C1NCBI Gene ID:63894
Calculated MW	57kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using VIPAS39 antibody at 1:1000 dilution.

Background

This gene encodes a protein involved in the sorting of lysosomal proteins. Mutations in this gene are associated with ARCS2 (arthrogryposis, renal dysfunction, and cholestasis-2). Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.