

ERO1LB Polyclonal Antibody

Catalog No: #30477

Package Size: #30477-1 50ul #30477-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ERO1LB Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms
Immunogen Description	Recombinant fusion protein of human ERO1LB (NP_063944.3).
Other Names	ERO1B; ERO1LB; Ero1beta; ERO1-like protein beta
Accession No.	Swiss-Prot#:Q86YB8NCBI Gene ID:56605
Calculated MW	48kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ERO1LB antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.