

OSBP Polyclonal Antibody

Catalog No: #30510



Package Size: #30510-1 50ul #30510-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	OSBP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human OSBP (NP_002547.1).
Other Names	OSBP;OSBP1
Accession No.	Uniprot:P22059GeneID:5007
Calculated MW	110kDa
SDS-PAGE MW	110kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using OSBP antibody.



Immunofluorescence analysis of NIH-3T3 cells using OSBP Rabbit pAb.

Background

Oxysterol binding protein is an intracellular protein that is believed to transport sterols from lysosomes to the nucleus where the sterol down-regulates the genes for the LDL receptor, HMG-CoA reductase, and HMG synthetase

Note: This product is for in vitro research use only and is not intended for use in humans or animals.