

## RSK1 p90 (Phospho-Ser380) Antibody

Catalog No: #12819



Package Size: #12819-1 50ul #12819-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RSK1 p90 (Phospho-Ser380) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purification using immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	RSK1 p90 (phospho S380) Antibody detects endogenous levels of RSK1 p90 only when phosphorylated at Ser380
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human RSK1 p90 around the phosphorylation site of Ser380.
Target Name	RSK1 p90 (Phospho-Ser380)
Modification	Phospho
Other Names	RPS6KA1, MAPKAPK-1a, p90RSK, p90-RSK 1, Pp90RSK1, Ribosomal S6 kinase 1, RSK-1, S6K-alpha 1, MAPKAP kinase 1a, MAPKAPK1A, p90 s6 kinase, p90S6K, RSK, RSK1, HU-1, p90RSK1, S6K-alpha-1
Accession No.	Swiss-Prot#: Q15418NCBI Gene ID: 6195
Target Species	human
Calculated MW	90 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.