

TUBGCP3 Polyclonal Antibody

Catalog No: #30561

Package Size: #30561-1 50ul #30561-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TUBGCP3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human TUBGCP3 (NP_001273207.1).
Other Names	TUBGCP3;104p;GCP3;Grip104;SPBC98;Spc98;Spc98p
Accession No.	Uniprot:Q96CW5GeneID:10426
Calculated MW	105kDa
SDS-PAGE MW	105kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TUBGCP3 antibody.

Background

Gamma-tubulin complex is necessary for microtubule nucleation at the centrosome.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.