

EXOSC8 Polyclonal Antibody

Catalog No: #30578

Package Size: #30578-1 50ul #30578-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	EXOSC8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu
Immunogen Description	Recombinant fusion protein of human EXOSC8 (NP_852480.1).
Other Names	EXOSC8; CIP3; EAP2; OIP2; PCH1C; RRP43; Rrp43p; bA421P11.3; p9; exosome component 8
Accession No.	Swiss-Prot#:Q96B26NCBI Gene ID:11340
Calculated MW	35kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using EXOSC8 antibody.

Background

This gene encodes a 3'-5' exoribonuclease that specifically interacts with mRNAs containing AU-rich elements. The encoded protein is part of the exosome complex that is important for the degradation of numerous RNA species. A pseudogene of this gene is found on chromosome 6.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.