

EIF5A2 Polyclonal Antibody

Catalog No: #30619



Package Size: #30619-1 50ul #30619-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	EIF5A2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human EIF5A2 (NP_065123.1).
Other Names	EIF5A2; EIF-5A2; eIF5AII; eukaryotic translation initiation factor 5A2
Accession No.	Swiss-Prot#:Q9GZV4NCBI Gene ID:56648
Calculated MW	17kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:200 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using EIF5A2 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.