

MYL7 Polyclonal Antibody

Catalog No: #30621

Package Size: #30621-1 50ul #30621-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MYL7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human MYL7 (NP_067046.1).
Other Names	MYL7; MYL2A; MYLC2A; myosin light chain 7
Accession No.	Swiss-Prot#:Q01449NCBI Gene ID:58498
Calculated MW	19kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using MYL7 antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.