

HNF-4 α/γ (Acetyl Lys127/79) Polyclonal Antibody

Catalog No: #HW179



Package Size: #HW179-1 50ul #HW179-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	HNF-4 α/γ (Acetyl Lys127/79) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous HNF-4 α/γ when Acetyl occurs at Lys127B α 79(gamma)
Immunogen Type	protein
Immunogen Description	Synthetic Acetyl peptide from human protein at AA range: 127(HNF-4 α)/79(HNF-4 γ)
Target Name	HNF-4 α/γ (Acetyl Lys127/79)
Other Names	Hepatocyte nuclear factor 4-alpha (HNF-4-alpha) (Nuclear receptor subfamily 2 group A member 1) (Transcription factor 14) Hepatocyte nuclear factor 4-gamma (HNF-4-gamma) (Nuclear receptor subfamily 2 group A member 2)(TCF-14) (Transcription factor HNF-4)
Accession No.	Swiss-Prot:P41235/Q14541NCBI Gene ID:3172
Calculated MW	55 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B α 1500-1B β 12000

Images



Western blot analysis of mouse-lung mouse-brain mouse-heart Hela mouse-liver lysate, antibody was diluted at 500. Secondary antibody was diluted at 1:20000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.