

MARCKSL1 Polyclonal Antibody

Catalog No: #30833

Package Size: #30833-1 50ul #30833-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MARCKSL1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human MARCKSL1 (NP_075385.1).
Other Names	MARCKSL1;F52;MACMARCKS;MLP;MLP1;MRP
Accession No.	Uniprot:P49006GeneID:65108
Calculated MW	35kDa
SDS-PAGE MW	35kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using MARCKSL1 antibody.

Background

This gene encodes a member of the myristoylated alanine-rich C-kinase substrate (MARCKS) family. Members of this family play a role in cytoskeletal regulation, protein kinase C signaling and calmodulin signaling. The encoded protein affects the formation of adherens junction. Alternative splicing results in multiple transcript variants. Pseudogenes of this gene are located on the long arm of chromosomes 6 and 10.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.