

ELMO3 Polyclonal Antibody

Catalog No: #30842

Package Size: #30842-1 50ul #30842-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ELMO3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ELMO3 (NP_078988.2).
Other Names	ELMO3;CED-12;CED12;ELMO-3
Accession No.	Uniprot:Q96BJ8GeneID:79767
Calculated MW	81kDa
SDS-PAGE MW	81kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ELMO3 antibody.

Background

The protein encoded by this gene is similar to a *C. elegans* protein that functions in phagocytosis of apoptotic cells and in cell migration. Other members of this small family of engulfment and cell motility (ELMO) proteins have been shown to interact with the dedicator of cyto-kinesis 1 protein to promote phagocytosis and effect cell shape changes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.