

## DPF1 Polyclonal Antibody

Catalog No: #30897



Package Size: #30897-1 50ul #30897-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	DPF1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human DPF1 (NP_001128627.1).
Other Names	DPF1; BAF45b; NEUD4; neuro-d4; double PHD fingers 1
Accession No.	Swiss-Prot#:Q92782NCBI Gene ID:8193
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:100

## Images



Immunofluorescence analysis of A549 cells using DPF1 antibody.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.