

ZNF544 Polyclonal Antibody

Catalog No: #30904



Package Size: #30904-1 50ul #30904-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ZNF544 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human ZNF544 (NP_055295.2).
Other Names	ZNF544; zinc finger protein 544
Accession No.	Swiss-Prot#:Q6NX49NCBI Gene ID:27300
Calculated MW	82kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZNF544 antibody at 1:1000 dilution.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.