

ZSCAN26 Polyclonal Antibody

Catalog No: #30968

Package Size: #30968-1 50ul #30968-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZSCAN26 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human ZSCAN26 (NP_689949.3).
Other Names	ZSCAN26;SRE-ZBP;SREZBP;ZNF187
Accession No.	Uniprot:Q16670GeneID:7741
Calculated MW	38kDa
SDS-PAGE MW	38kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse brain, using ZSCAN26 antibody.

Background

May be involved in transcriptional regulation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.