

TAMM41 Polyclonal Antibody

Catalog No: #31530

Package Size: #31530-1 50ul #31530-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | TAMM41 Polyclonal Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Affinity purification |
| Applications | WB |
| Species Reactivity | Human,Mouse |
| Immunogen Description | Recombinant fusion protein of human TAMM41 (NP_620162.1). |
| Other Names | TAMM41;C3orf31;RAM41;TAM41 |
| Accession No. | Uniprot:Q96BW9GeneID:132001 |
| Calculated MW | 38kDa |
| SDS-PAGE MW | 38kDa |
| Formulation | PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| Storage | Store at -20°C. Avoid freeze / thaw cycles. |

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TAMM41 antibody.

Background

Catalyzes the conversion of phosphatidic acid (PA) to CDP-diacylglycerol (CDP-DAG, an essential intermediate in the synthesis of phosphatidylglycerol, cardiolipin and phosphatidylinositol).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.