

## TGFB111 Polyclonal Antibody

Catalog No: #31555

Package Size: #31555-1 50ul #31555-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	TGFB111 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human TGFB111 (NP_001035919.1).
Other Names	TGFB111;ARA55;HIC-5;HIC5;TSC-5
Accession No.	Uniprot:O43294GeneID:7041
Calculated MW	50kDa
SDS-PAGE MW	50KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of Mouse lung, using TGFB111 antibody.



Western blot analysis of extracts of various cell lines, using TGFB111 antibody.

## Background

This gene encodes a coactivator of the androgen receptor, a transcription factor which is activated by androgen and has a key role in male sexual differentiation. The encoded protein is thought to regulate androgen receptor activity and may have a role to play in the treatment of prostate cancer. Multiple transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.