

VEPH1 Polyclonal Antibody

Catalog No: #31581

Package Size: #31581-1 50ul #31581-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	VEPH1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human VEPH1 (NP_078897.2).
Other Names	VEPH1; MELT; VEPH; ventricular zone expressed PH domain containing 1
Accession No.	Swiss-Prot#:Q14D04NCBI Gene ID:79674
Calculated MW	94kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using VEPH1 antibody at 1:1000 dilution.



Western blot analysis of extracts of Mouse brain, using VEPH1 antibody at 1:1000 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.