

MTERFD3 Polyclonal Antibody

Catalog No: #31585

Package Size: #31585-1 50ul #31585-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MTERFD3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human MTERFD3 (NP_079474.2).
Other Names	MTERF2;MTERFD3;mTERFL
Accession No.	Uniprot:Q49AM1GeneID:80298
Calculated MW	38-44kDa
SDS-PAGE MW	44KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using MTERFD3 antibody.

Background

Binds mitochondrial DNA and plays a role in the regulation of transcription of mitochondrial mRNA and rRNA species.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.