

## EXOSC4 Polyclonal Antibody

Catalog No: #31675

Package Size: #31675-1 50ul #31675-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	EXOSC4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Hu,Rt
Immunogen Description	Recombinant fusion protein of human EXOSC4 (NP_061910.1).
Other Names	EXOSC4; RRP41; RRP41A; Rrp41p; SKI6; Ski6p; hRrp41p; p12A; exosome component 4
Accession No.	Swiss-Prot#:Q9NPD3NCBI Gene ID:54512
Calculated MW	26kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000 IHC □ 1:100 - 1:200

## Images



Western blot analysis of extracts of various cell lines, using EXOSC4 antibody.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.