

EBNA1BP2 Polyclonal Antibody

Catalog No: #31676



Package Size: #31676-1 50ul #31676-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	EBNA1BP2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human EBNA1BP2 (NP_006815.2).
Other Names	EBNA1BP2; EBP2; NOBP; P40; EBNA1 binding protein 2
Accession No.	Swiss-Prot#:Q99848NCBI Gene ID:10969
Calculated MW	40kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:1000 - 1:2000 IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using EBNA1BP2 antibody at 1:1000 dilution.



Confocal immunofluorescence analysis of U2OS cells using EBNA1BP2 Polyclonal antibody at dilution of 1:200. Blue: DAPI for nuclear staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.