

## NFIX Polyclonal Antibody

Catalog No: #31707

Package Size: #31707-1 50ul #31707-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	NFIX Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms
Immunogen Description	Recombinant fusion protein of human NFIX (NP_002492.2).
Other Names	NFIX; CTF; MRSHSS; NF-I/X; NF1-X; NF1A; SOTOS2; nuclear factor I X
Accession No.	Swiss-Prot#:Q14938NCBI Gene ID:4784
Calculated MW	50kDa
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:1000

## Images



Western blot analysis of extracts of various cell lines, using NFIX antibody.

## Background

The protein encoded by this gene is a transcription factor that binds the palindromic sequence 5'-TTGGCNNNNNGCCAA-3 in viral and cellular promoters. The encoded protein can also stimulate adenovirus replication in vitro. Three transcript variants encoding different isoforms have been found for this gene.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.