

ZRANB3 Polyclonal Antibody

Catalog No: #31722

Package Size: #31722-1 50ul #31722-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ZRANB3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human ZRANB3 (NP_001273497.1).
Other Names	ZRANB3; 4933425L19Rik; AH2; zinc finger RANBP2-type containing 3
Accession No.	Swiss-Prot#:Q5FWF4NCBI Gene ID:84083
Calculated MW	123kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:1000 - 1:3000

Images



Western blot analysis of extracts of various cell lines, using ZRANB3 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.