

EXOSC7 Polyclonal Antibody

Catalog No: #31733



Package Size: #31733-1 50ul #31733-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | EXOSC7 Polyclonal Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Affinity purification |
| Applications | WB |
| Species Reactivity | Hu,Ms,Rt |
| Immunogen Description | Recombinant fusion protein of human EXOSC7 (NP_055819.2). |
| Other Names | EXOSC7; EAP1; RRP42; Rrp42p; hRrp42p; p8; exosome component 7 |
| Accession No. | Swiss-Prot#:Q15024NCBI Gene ID:23016 |
| Calculated MW | 37kDa |
| Formulation | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4. |
| Storage | Store at -20°C |

Application Details

WB 1:200 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using EXOSC7 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.