

# YWHAZ Conjugated Antibody

Catalog No: #C28706

Package Size: #C28706-AF350 100ul #C28706-AF405 100ul #C28706-AF488 100ul #C28706-AF555 100ul #C28706-AF594 100ul #C28706-AF647 100ul #C28706-AF680 100ul #C28706-AF750 100ul #C28706-Biotin 100ul #C28706-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	YWHAZ Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB, IF
Species Reactivity	Hu,Ms,Rt
Immunogen Description	A synthetic peptide of human YWHAZ (NP_001129171.1).
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	YWHAZ; 14-3-3-zeta; HEL-S-3; HEL-S-93; HEL4; KCIP-1; YWHAD; 14-3-3 protein zeta/delta
Accession No.	Swiss-Prot#:P63104NCBI Gene ID:7534
Calculated MW	Refer to figures
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Background

This gene product belongs to the 14-3-3 family of proteins which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 99% identical to the mouse, rat and sheep orthologs. The encoded protein interacts with IRS1 protein, suggesting a role in regulating insulin sensitivity. Several transcript variants that differ in the 5' UTR but that encode the same protein have been identified for this gene.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.