

TP53RK Conjugated Antibody

Catalog No: #C28774



Package Size: #C28774-AF350 100ul #C28774-AF405 100ul #C28774-AF488 100ul #C28774-AF555 100ul #C28774-AF594 100ul #C28774-AF647 100ul #C28774-AF680 100ul #C28774-AF750 100ul #C28774-Biotin 100ul #C28774-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TP53RK Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB, IF
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human TP53RK (NP_291028.3).
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	TP53RK; BUD32; C20orf64; Nori-2; Nori-2p; PRPK; dJ101A2; TP53-regulating kinase
Accession No.	Swiss-Prot#:Q96S44NCBI Gene ID:112858
Calculated MW	28kDa
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.