

## MYH6 Polyclonal Antibody

Catalog No: #46918



Package Size: #46918-1 50ul #46918-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	MYH6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MYH6 Polyclonal Antibody detects endogenous levels of MYH6
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human MYH6. at AA range: 341-390
Target Name	MYH6
Other Names	Myosin-6 (Myosin heavy chain 6) (Myosin heavy chain, cardiac muscle alpha isoform) (MyHC-alpha)
Accession No.	Swiss-Prot:P13533NCBI Gene ID:4624
Calculated MW	200 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1B£B'1500-1B£B'12000

## Images



Western Blot analysis of mouse-heart cells using primary antibody diluted at 1:2000(°C overnight). Secondary antibodyB£B'Goat Anti-rabbit IgG IRDye 800( diluted at 1:5000, 25B'£C, 1 hour)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.