

YTHDF3 Conjugated Antibody

Catalog No: #C31538

Package Size: #C31538-AF350 100ul #C31538-AF405 100ul #C31538-AF488 100ul #C31538-AF555 100ul #C31538-AF594 100ul #C31538-AF647 100ul #C31538-AF680 100ul #C31538-AF750 100ul #C31538-Biotin 100ul #C31538-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	YTHDF3 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB, IF
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human YTHDF3 (NP_689971.4).
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	YTHDF3; YTH domain-containing family protein 3
Accession No.	Swiss-Prot#:Q7Z739NCBI Gene ID:253943
Calculated MW	64kDa
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Background

This gene encodes a member of the YTH (YT521-B homology) domain protein family. The YTH domain is common in eukaryotes, is often found in the middle of the protein sequence, and may function in binding to RNA. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.