

# ZNF483 antibody

Catalog No: #31828



Package Size: #31828-1 50ul #31828-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZNF483 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Immunogen Description	Full length fusion protein
Target Name	ZNF483
Other Names	ZSCAN48; ZKSCAN16
Accession No.	NCBI Protein#:BC065738
Concentration	0.5mg/ml
Formulation	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Storage	Store at -20°C/1 year

## Application Details

IHC dilution:1:25-1:100

## Images



The image is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using (ZNF483 Antibody) at dilution 1/30.

## Background

May be involved in transcriptional regulation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.