

## DHRS7B antibody

Catalog No: #31902



Package Size: #31902-1 50ul #31902-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	DHRS7B antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Immunogen Description	Fusion protein of human DHRS7B
Target Name	DHRS7B
Other Names	CGI-93; SDR32C1
Accession No.	NCBI Protein#:BC004126
Concentration	1.4mg/ml
Formulation	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Storage	Store at -20°C/1 year

## Application Details

IHC dilution:1:50-1:100

## Images



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using (Antibody) at dilution 1/80.

## Background

This gene is located within the Smith-Magenis syndrome region on chromosome 17. It encodes a protein of unknown function.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.