

NHLRC2 antibody

Catalog No: #31912

Package Size: #31912-1 50ul #31912-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NHLRC2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Immunogen Description	Fusion protein of human NHLRC2
Target Name	NHLRC2
Accession No.	NCBI Protein#:BC032598
Concentration	0.6mg/ml
Formulation	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at -20°C/1 year

Application Details

IHC dilution:1:50-1:200

Images



The image is immunohistochemistry of paraffin-embedded Human breast cancer tissue using (NHLRC2 Antibody) at dilution 1/40.

Background

The specific function of this protein remains unknown.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.