

ZDHHC3 antibody

Catalog No: #31998

Package Size: #31998-1 50ul #31998-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZDHHC3 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu, Ms
Immunogen Description	Fusion protein of human ZDHHC3
Target Name	ZDHHC3
Other Names	GODZ; DHHC3; DHHC-3; ZNF373
Accession No.	NCBI Protein#:BC015467
Concentration	1mg/ml
Formulation	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at -20°C/1 year

Application Details

IHC dilution:1:50-1:100

Images



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using (Zntibody) at dilution 1/35.

Background

Palmitoyltransferase with broad specificity (PubMed:21926431, PubMed:23034182). Palmitoylates GABA receptors on their gamma subunit (GABRG1, GABRG2 and GABRG3), which regulates synaptic clustering and/or cell surface stability (By similarity). Palmitoylates glutamate receptors GRIA1 and GRIA2, which leads to their retention in Golgi. May also palmitoylate DLG4, DNAJC5 and SNAP25 (By similarity).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.