

## DIAPH2 antibody

Catalog No: #43485

Package Size: #43485-1 50ul #43485-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	DIAPH2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu, Ms
Immunogen Description	Fusion protein of human DIAPH2
Target Name	DIAPH2
Other Names	DIA; POF; DIA2; DRF2; POF2; POF2A
Accession No.	NCBI Protein#:BC117414
Concentration	0.8mg/ml
Formulation	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Storage	Store at -20°C/1 year

## Application Details

IHC dilution:1:50-1:100

## Images



The image is immunohistochemistry of paraffin-embedded Human breast cancer tissue using (Antibody) at dilution 1/55.

## Background

The product of this gene belongs to the diaphanous subfamily of the formin homology family of proteins. This gene may play a role in the development and normal function of the ovaries. Defects in this gene have been linked to premature ovarian failure 2. Alternatively spliced transcript variants encoding different isoforms have been identified.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.