

ZNF35 antibody

Catalog No: #43490

Package Size: #43490-1 50ul #43490-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF35 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Immunogen Description	Synthetic peptide of human ZNF35
Target Name	ZNF35
Other Names	HF10; HF.10; Zfp105
Accession No.	NCBI Protein#:NP_003411
Concentration	1.1mg/ml
Formulation	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at -20°C/1 year

Application Details

IHC dilution:1:50-1:100

Images



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using (ZNF35 Antibody) at dilution 1/40.

Background

May be involved in transcriptional regulation. Involved in cell differentiation and/or proliferation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.