

DPCD antibody

Catalog No: #43496

Package Size: #43496-1 50ul #43496-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DPCD antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu, Ms, Rt
Immunogen Description	Synthetic peptide of human DPCD
Target Name	DPCD
Accession No.	NCBI Protein#:NP_056263
Concentration	1.5mg/ml
Formulation	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at -20°C/1 year

Application Details

IHC dilution:1:50-1:100

Images



The image is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using (tibody) at dilution 1/50.

Background

This gene in mouse encodes a protein that may be involved in the generation and maintenance of ciliated cells. In mouse, expression of this gene increases during ciliated cell differentiation, and disruption of this gene has been linked to primary ciliary dyskinesia.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.