

## HSDL2 Conjugated Antibody

Catalog No: #C31905



Package Size: #C31905-AF350 100ul #C31905-AF405 100ul #C31905-AF488 100ul #C31905-AF555 100ul #C31905-AF594 100ul  
 #C31905-AF647 100ul #C31905-AF680 100ul #C31905-AF750 100ul #C31905-Biotin 100ul #C31905-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	HSDL2 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB, IF
Species Reactivity	Hu
Immunogen Description	Fusion protein of human HSDL2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	HSDL2
Other Names	C9orf99; SDR13C1
Accession No.	Swiss-Prot#: Q6YN16NCBI Protein#: BC004331
Calculated MW	45 kDa
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Background

Has apparently no steroid dehydrogenase activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.