

Rat HB- α 1 ELISA Kit

Catalog No: #EK20109



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat HB- α 1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rt

Application Details

Detect Range: 0.1 - 16 ng/mL

Sample Type: Serum, plasma, Urine, tissue homogenates and other biological fluids.

Sample Volume: 50 μ l

Assay Time: 3-5h

Detection wavelength: 450nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.