

# Human Melanocortin receptor 4 (MC4R) ELISA Kit

Catalog No: #EK20115



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                   |  |
|-------------------|--|
| Product Name      | Human Melanocortin receptor 4 (MC4R) ELISA Kit |
| Brief Description | Elisa Kit                                      |
| Applications      | ELISA  |
| Specificity       | Homo sapiens human                             |
| Other Names       | Melanocortin receptor 4, MC4-R, MC4R           |
| Accession No.     | UniProtKB:P32245 Genes:160                     |
| Storage           | 12 months                                      |

## Application Details

|                       |  |
|-----------------------|--|
| Detect Range:         | 0.156-10ng/mL  |
| Sensitivity:          | 0.060ng/mL   |
| Sample Type:          | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Sample Volume:        | 100 µL   |
| Assay Time:           | 1-4.5h   |
| Detection wavelength: | 450 nm   |

## Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.