

Recombinant Human Novel Coronavirus Spike glycoprotein(S), partial

Catalog No: #AP89528

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Novel Coronavirus Spike glycoprotein(S), partial
Brief Description	Recombinant Protein
Purification	Greater than 85% as determined by SDS-PAGE.
Other Names	Spike glycoprotein
Accession No.	Uniprot ID: MN908947
Calculated MW	78.0kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	VNLTRTQLPPAYTNSFTRGVYYPDKVFRSSVLHSTQDLFLPFFSNVTWFHAIHVSGTNGTKRFDNPVLPFND GVYFASTSEKSNIRGWIFGTTLDSKTQSLIVNATNVVIK/CEFQFCNDPFLGVYHKNKSWMESEFRVYSS ANNCTFEYVSQPFLMDLEGKQGNFKNLREFVFKNIDGYFKIYSKHTPINLVRDLPQGFSALEPLVDLPIGINITRF QTLALHRSYLTPGDSSSGWTAGAAAYVGYLQPRFTLLKYNENGTITDAVDCALDPLSETKCTLSFTVEKGI YQTSNFRVQPTESIVRFPNITNLCPFGEVFNATRFASVYAWNKRKISNCVADYSVLVNSASFSTFKCYGVSPK LNDLCFTNVYADSFVIRGDEVQRQIAPGQTGKIADYNYKLPDFTGCVIAWNSNLDKSKVGGNYNYLRLFRKS NLKPFERDISTEIYQAGSTPCNGVEGFNCYFPLQSYGFQPTNGVGYQPYRVVLSFELLHAPATVCGPKKSTN LVKNKCVNFNGLTGTGVLTESNKKFLPFQFGQGRDIADTTDAVRDPQTLEILDITPCSFGGVSVITPGTNTSN QVAVLYQDVNCTEVPVAIHADQLTPTWRVYSTGSNVFQTRAGCLIGAEHVNSYECDIPIGAGICASYQTQTN SPRRAR
Storage	Store at -20°C for extended storage, conserve at -20°C or -80°C

Application Details

Expression Region: 16-685aa

Sequence Info: partial Spike protein S1 chain

Research Topic: Microbiology

Organism: Human Novel Coronavirus (SARS-CoV-2/ 2019-nCoV)

Storage Buffer:

Product Description

Repeated freezing and thawing is not recommended. Store working aliquots at -20°C for up to one week.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.