

Goat Anti-Chicken IgY (H&L) HRP

Catalog No: #L35001

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Chicken IgY (H&L) HRP
Brief Description	Conjugated Secondary Antibody
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH, IC
Species Reactivity	C
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Chicken IgY. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Chicken IgY
Conjugates	HRP
Concentration	1 mg/ml
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

Published Papers

el at., SynCAM1 deficiency in the hippocampal parvalbumin interneurons contributes to sevoflurane-induced cognitive impairment in neonatal rats
CNS Neurosci Ther On2023 Dec 17byMing-Ming Zhao , Han-Wen Gu et al..PMID:38105652, , (2023)

[PMID:38105652](https://pubmed.ncbi.nlm.nih.gov/38105652/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.