

Goat Anti-Rabbit IgG (H&L) HRP

Catalog No: #L35009



Package Size: #L35009-1 100 ul #L35009-2 500 ul #L35009-3 1 ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Rabbit IgG (H&L) HRP
Brief Description	Conjugated Secondary Antibody
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH, IC
Species Reactivity	Rb
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Rabbit IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Rabbit IgG
Conjugates	HRP
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

Published Papers

el at., TSC22D3 as an immune-related prognostic biomarker for acute myeloid leukemia. *Science* 2023 Jul 22;381(6652):1215-1220. doi: 10.1126/science.1258112. PMID: 37575189, (2023)

[PMID:37575189](https://pubmed.ncbi.nlm.nih.gov/37575189/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.