

# Rabbit Anti-Pig IgG (H&L) HRP

Catalog No: #L35026



Package Size: #L35026-1 100 ul #L35026-2 500 ul #L35026-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |  |
|-----------------------|--|
| Product Name          | Rabbit Anti-Pig IgG (H&L) HRP  |
| Brief Description     | Conjugated Secondary Antibody  |
| Host Species          | Rabbit   |
| Clonality             | Polyclonal   |
| Purification          | The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.                                       |
| Applications          | E, WB, IH, IC  |
| Species Reactivity    | P  |
| Specificity           | By immunoelectrophoresis and ELISA this antibody reacts specifically with Pig IgG. No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Pig IgG  |
| Conjugates            | HRP  |
| Formulation           | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide   |

## Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.