

## Rabbit Anti-Rat IgG (H&amp;L) HRP

Catalog No: #L35027



Package Size: #L35027-1 100 ul #L35027-2 500 ul #L35027-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rabbit Anti-Rat IgG (H&L) HRP
Brief Description	Conjugated Secondary Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH, IC
Species Reactivity	R
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Rat IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Rat IgG
Conjugates	HRP
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

## Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

## Published Papers

el at., Klotho in Osx+-mesenchymal progenitors exerts pro-osteogenic and anti-inflammatory effects during mandibular alveolar bone formation and repair. In Signal Transduct Target Ther on 2022 May 11 by Yi Fan, Chen Cui, et al..PMID: 35538062, , (2022)

[PMID:35538062](https://pubmed.ncbi.nlm.nih.gov/35538062/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.