

# Mouse Anti-Chicken IgY (H&L) HRP

Catalog No: #L35029



Package Size: #L35029-1 100 ul #L35029-2 500 ul #L35029-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Anti-Chicken IgY (H&L) HRP
Brief Description	Conjugated Secondary Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH, IC
Species Reactivity	C
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Chicken IgY. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Chicken IgY
Conjugates	HRP
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

## Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.