

# Rabbit Anti-Monkey IgG (H&L) Biotin

Catalog No: #L35056



Package Size: #L35056-1 100 ul #L35056-2 500 ul #L35056-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rabbit Anti-Monkey IgG (H&L) Biotin
Brief Description	Conjugated Secondary Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH, IC
Species Reactivity	Mk
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Monkey IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Monkey IgG
Conjugates	Biotin
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

## Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.