

Rabbit Anti-Mouse IgG (H&L) Biotin

Catalog No: #L35057



Package Size: #L35057-1 100 ul #L35057-2 500 ul #L35057-3 1 ml

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rabbit Anti-Mouse IgG (H&L) Biotin
Brief Description	Conjugated Secondary Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH, IC
Species Reactivity	M
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Mouse IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Mouse IgG
Conjugates	Biotin
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.