

# Rabbit Anti-Sheep IgG (H&L) Biotin

Catalog No: #L35060



Package Size: #L35060-1 100 ul #L35060-2 500 ul #L35060-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rabbit Anti-Sheep IgG (H&L) Biotin
Brief Description	Conjugated Secondary Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH, IC
Species Reactivity	S
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Sheep IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Sheep IgG
Conjugates	Biotin
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

## Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

## Published Papers

el at., The relationship between myodural bridge, atrophy and hyperplasia of the suboccipital musculature, and cerebrospinal fluid dynamicsInSci RepOn2023 Nov 2byHeng Yang # 1, Xiao-Song Wei et al..PMID: 37919345, , (2023)  
[PMID:37919345](https://pubmed.ncbi.nlm.nih.gov/37919345/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.